

Fingerprinting : the new frontier – breaking barriers to implementation

Prof. **Simona Francese**

CMSI, BMRC, Sheffield Hallam University, Sheffield, UK

s.francese@shu.ac.uk



Fingerprinting



Yorkshire and Humber Regional Scientific Support Services

Courtesy of Peter Arnold Head of Operations YHRSSS




Conventional fingerprinting pitfalls

The "soul" of a fingerprint



Chemical criminal profiling

Francesse S, et al. *Analyst*. 2013, 142(14):2518-46
Francesse S et al. *Analyst* 2017, 138(15): 4215-28



Chemical criminal profiling

Mass Spectrometry (and hyphenated techniques)

DESI
DART
MALDI
Ag LDI
Paper Spray
SALDI
SIMS

LC ESI MS
GC ESI MS




Spectroscopy

Raman
FT IR and ATR FT IR


Imaging and profiling capabilities

Number of publications on fingerprint analysis per analytical technique (2008-2018)




Francese S, 2019, Austr. J. Forensic Sci. 51 (6): 623-635

WHAT IS MASS SPECTROMETRY (MS)?



Scales...for molecules!!!!



Mass Spectrometry as investigative analytical tool



Before Christmas



After Christmas



caffeine weighs =
194.1 Daltons




Before Christmas

After Christmas




caffeine weighs =
194.1 Daltons


MALDI MS Profiling OF LATENT FINGERMARKS



MALDI MSP and MALDI MSI OF LATENT FINGERMARKS




MALDI MSI OF LATENT FINGERMARKS




PHYSICAL INFORMATION

Improving the coverage of the ridge pattern







CHEMICAL INFORMATION


(very few examples)

Personal information (*endogenous molecules*)

CAN WE DETECT SEX FROM FINGERPRINTS?




Males and Females are discriminated with the **75% of confidence** (85% if less harsh classification criteria were applied)



Personal information (*endogenous molecules*)

XGBOOST classification with full consensus scoring

86% accuracy of prediction – considered by the Police viable for triaging crime scene marks




Application of “full consensus scoring” strategy. Correlation between calculated majority vote scores (x-axis) and full consensus scoring schemes (y-axis), for the XGBOOST classifier.

Each scatter point corresponds to the average cross-validation accuracy score for a specific peak picking strategy, and k -fold, with each of the four feature selection strategies treated as separate scatter points. The diagonal line $y = x$ is shown.


Pharmacological states (*semi-endogenous molecules*)

Propranolol

m/z 260.1650




- **4- hydroxypropranolol**




m/z 276.1599

- **N- desisopropylpropranolol**



m/z 217.1103

Confirmation of presence of 4-hydroxypropranolol in fingerprint through drift time assessment of its product ion (m/z 199.0940)




Assessing crime dynamics (and not only) (*exogenous molecules*)

**MALDI MSP/MSI
combined proteomic
approach for the
specific, reliable and
informative analysis
of blood evidence**



9 YEAR OLD AB-1 treated BLOOD PALM PRINT



PRINT




Patel et al, 2015, *Analyst*, 141(1):191-8

Francesse S. 2019, *Austr. J. Forensic Sci.* 51 (6): 623-635

37 YEAR OLD NIN treated blood mark on fabric




Deininger et al, 2016, *Proteomics* 16(11-12):1707-17




MALDI MSI of animal blood signatures

4 YEAR OLD AB-1 treated bovine blood mark



Particularly the images of the ions at m/z 1669.832, 784.434 and 1477.798 depict the partial ridge detail that was observed in the optical images (red rectangle)

MALDI MSI of overlapping pre-enhanced fingerprints contaminated with HbC and HbJ-Baltimore blood variants



Breaking barriers to implementation

Understanding the use and limitations of the technology

Understanding end user needs

Disseminate/increase awareness

Walking in the right direction to break barriers to implementation

Make research borderless

Actively involve end users in the developments

Breaking barriers to implementation

Understanding the use and limitations of the technology

Understanding end user needs

Disseminate/increase awareness

Walking in the right direction to break barriers to implementation


Make research borderless

Actively involve end users in the developments

To use MALDI....the fingerprint needs to be visualised first....Duh!!!




MALDI matrix



The dry-wet method

Patent no. GB2489215; EU application (EP268401A1); Granted HK patent (HK1169854A1); Granted US patent (US9261438B)


LS Ferguson, R Wolstenholme, S. Francese



mark enhancement by matrix powdering



UV image



fluorescent microscope image




MALDI MSI image

The dry wet method adaptation – ENFSI collaborative exercise 2020




Compostable bag




Illumination at 365 nm




Powdering with MALDI Matrix



Illumination at 365 nm




MALDI Matrix solution spray



m/z 460.8


MALDI MS image

Breaking barriers to implementation



CSI place lifts on acetate for transport...would MALDI still work when the lift is peeled off???


Deposition of fingerprint onto aluminium plate




Aluminium powder



Primary lift



Primary lift peeled from acetate



Triacylglyceride
 m/z 664.6

Breaking barriers to implementation

Understanding the use and limitations of the technology


Understanding end user needs

Disseminate/increase awareness


Walking in the right direction to break barriers to implementation

Make research borderless

Actively involve end users in the developments



Breaking barriers to implementation




April 2014




Aug 2011

Fingerprint breakthrough offers new forensic evidence

A technology to extract fingerprints from a crime scene could show if a criminal suspect has taken drugs or been in contact with explosives.
 Researchers at Sheffield Hallam University have developed a way to capture fingerprints that can identify substances touched, and secretions from a person's body.
 This could provide new information about a suspect's actions and habits.



Mar 2014



BBC Click reports on the technology that is helping police learn more about people from their fingerprints.

Aug 2021

Home » Blog » 2021 08
 At the frontier of fingermarks:
 Prof Simona Francese on 13
 years of research into
 molecular fingerprinting



1. Fingermarks match a mark left at a crime scene - molecular fingerprinting
 2. Identifying if the person touched products or other materials, to blood.
 Professor Simona Francese, an expert in molecular and bioanalytical mass spectrometry, explains how she has spent 13 years researching how to identify fingermarks left at a crime scene for identification. Investigators take an impression of a fingermark, or wider impressions of a finger, and analyse the molecules - a fingerprint, so why is a fingerprint lies in its

ANZFSS 2018 [+ Follow](#)
 ANZFSS18 - Conference - 13th Sept 2018

Sept 2018



SIMONA FRANCESE
 Sheffield Hallam University
 BBC BREAKFAST **Oct 2017**




June 2014



June 2018

Breaking barriers to implementation....




Home Office funding 2011-2014 *“Integration of MALDI MSI into the Home Office operational fingerprint examination workflow”*





2013 - Attendance to crime scenes with West Yorkshire Police



- Understanding CSI requirements for collection of marks
- Understanding optimal storing of marks
- **Developing compatibility between fingerprint visualisation processing and MALDI**



Bradshaw et al FSI 2013 232(1-3):111-24



What Police thought in 2015.....

March
2015

Neil Denison

Head of identification services
and imaging services



'This technology shows real potential to be the next step change in developing fingerprint technology for the 21st century.'

Neil Denison

West Yorkshire Police

"This new fingerprinting technology has the potential to be a real step change in forensic work. This technology will allow us to gain significantly enhanced information from prints to build a picture of an offender's prior movements, and is now being trialled at crime scenes"

(Yorkshire Post **Aug 2015**)

Assistant Chief Constable, Russ Foster

West Yorkshire Police



Home Office Innovation Grant co-funded by West Yorkshire Police, UK 2014-2016
“Chemical Profiling of fingermarks : MALDI MSP and MALDI MSI method refinement for operational deployment in casework”


Analytical Methods

PAPER

Development of operational protocols for the analysis of primary and secondary fingermark lifts by MALDI-MS imaging†

Robert Bradshaw,^a Neil Denison^b and Simona Francese^{*a}

Eight years of intensive research have demonstrated that Matrix Assisted Laser Desorption Ionisation-Mass Spectrometry Profiling and Imaging (MALDI-MSP and MSI) are powerful tools to gather intelligence around a suspect lifestyle, directly from the identifying ridges of a latent fingermark. In the past three years, many efforts have been invested into translating laboratory methodologies to the field; this was undertaken by devising protocols either for (a) enabling initial fingermark visualisation, such as through the Dry-Wet method, recovery and subsequent MALDI MS based analysis, or for (b) rendering the MS methodologies




John Dixon, WYP


Operational casework (UK 2015)- harassment crime scene




Carbon Black
Lifted mark sealed in acetate




Cocaine In-Source Fragment
m/z 182.1



Cocaine
m/z 304.1



Home Office Recognition of the technology



Currently included
as a Cat C process
in the Fingermark
Visualisation
Manual

MALDI-MSI

Alternative Names
Matrix assisted laser desorption/ionisation-mass spectral imaging

Key Information

Where this process could be used
MALDI-MSI may be of use in identifying and mapping the distribution of unusual and/or significant contaminants present in fingerprints (e.g. drugs, explosives, condom lubricants etc.). The process may also be of use for 'filling in' ridge detail in marks that are only partially developed by conventional processes. See **Category B-C process options**.

Why the process is not in Category A
The full range of applications for MALDI-MSI has not yet been explored, nor has its compatibility with conventional development processes been explored in full.

Safety and Effectiveness Summary

The Process

- MALDI-MSI can only be carried out in a laboratory with specialist equipment and operator, which may have to be sourced within academic institutions.
- The effectiveness is linked to the compatibility of the fingerprint components and the matrix composition.

The Item or Surface

- The process is most effective at visualising latent marks on porous, semi-porous and adhesive surfaces although it can be used on non-porous surfaces.
- The process is most likely to require the mark to be lifted from the item as the chamber within the equipment is of very limited size. The lift is then processed.

Integrated Use
MALDI-MSI may be detrimental to subsequent fingerprint or forensic processing.

Further Reading

1. Bradshaw, R., Wolstenholme, R., Blackledge, R., Clench, M. R., Ferguson, L., Francese, S., Rapid Commun. Mass Spectrom. 2011, 3: 415-422.
2. Ferguson, L., Bradshaw, R., Wolstenholme, R., Clench, M. R., Francese, S., Anal. Chem. 2011, 83(14):5685-91

Competent personnel specialising in fingerprint visualisation must be consulted if considering the use of this process.

- Ensure all Category A process options have been explored before using this process.
- In the first instance, it is recommended that all relevant information is gathered from appropriate sources prior to considering the use of this process.

Process Overview
Matrix assisted laser desorption/ionisation-mass spectral imaging (MALDI-MSI) requires a homogeneous coating, known as a 'matrix', to first be applied to the fingerprint. Ions of a component present in latent and/or developed fingerprints are then desorbed from the surface by irradiation with a laser. The laser is moved between each firing so that the distribution of the ions can be mapped, giving an image of the fingerprint and/or chemical composition information.

**Cat C
TLR 3-4**

© See Photo Credits

promoted to Cat B
in the next DSTL
edition

Assessing Technology Readiness Levels

TRL	Description
1	Published papers reporting a technology or data potentially relevant to the area being considered
2	Published papers reporting a technology or data relevant to the area of application being considered
3	Initial application of technology being assessed to samples of forensic relevance in experimental settings
4	Focused studies showing feasibility of technology for the proposed application
5	Focused studies, assessing factors such as sensitivity and repeatability for the proposed application
6	Critical mass studies showing clear potential for operational use
7	Pseudo-operational trials, analysis of evidence under pseudo-operational conditions.
8	Operational trials on evidence obtained from casework
9	Publication of data from live casework
10	Publication of data showing the application of the technology in casework over several years

**TRL scale combined from FVM, Dstl,
EPSRC**

Assessing Technology Readiness Levels

Application of MALDI MSI to fingerprints	TRL	TRL description	Application of MALDI MS Profiling to fingerprints	TRL	TRL description	Recommendation
Visualisation of the ridge detail	9	Publication of data from live casework	NA	NA	NA	High Priority
Mapping of blood	7	Pseudo-operational trials, analysis of evidence under pseudo-operational conditions.	Detection and provenance of blood	8	Operational trials on evidence obtained from casework	High Priority
Mapping of drugs and metabolites	8	Operational trials on evidence obtained from casework	Detection and Identification of drugs and metabolites	8	Operational trials on evidence obtained from casework	High Priority
NA	NA	NA	Determination of sex from fingerprints	6	Critical mass studies showing clear potential for operational use	Medium priority
Mapping of condom lubricants	5	Focussed studies, assessing factors such as sensitivity and repeatability for the proposed application	Detection and identification of condom lubricants	8	Operational trials on evidence obtained from casework	Low priority
Mapping of pharmaceuticals and metabolites	2	Published papers reporting a technology or data relevant to the area of application being considered	Detection and identification of pharmaceuticals and metabolites	2	Published papers reporting a technology or data relevant to the area of application being considered	Low priority
Mapping of explosives	2	Published papers reporting a technology or data relevant to the area of application being considered	Detection and Identification of explosives	2	Published papers reporting a technology or data relevant to the area of application being considered	Low priority

So why isn't the technology used more often in major crimes?




1. Many Forensic providers and Major crime Managers lack of awareness of capabilities (and forgetfulness) despite major dissemination
 - No other way than keeping at it as much as possible!!
 - Round table with policing and new technologies/products representative discussing (non live) casework would be helpful (CAST used to organise this)
2. Live cases and backlogs don't lend themselves to decisions to explore new approaches
 - Tough one. Structured Programmes of CPD would help. This is also raised in the HoL report
3. Some confusion on exploitation
 - Are validation (and accreditation) obstacles to implementation? NOT AT AN INVESTIGATIVE LEVEL. Your thoughts in a court environment?

So why isn't the technology used more often in major crimes?

4. Innovation improves turnaround and employee satisfaction but may negate benefits of emerging capabilities
 - NO quick fix. Keep raising the issue at a management level. IN MAJOR CRIMES, If lifts are taken at the time, opportunity for additional intelligence will be maximised
5. Scepticism towards an emerging approach and fear to “burn” evidence.
 - Continue to offer evidence of operationally feasible technology

Pre-validation of MALDI MSP for the detection of blood and its provenance

Example 1




- **Access to at least hundreds of crime scene marks (non live cases) to establish the impact of :**

- Fingerprint molecular composition variability
- Fingerprint age
- Environmental factors
- Fingerprint visualisation processes
- Cross talk of the above factors

ON THE:

Success rate of this technology for the provision of profiling and biometric information

Research Councils...these mythological creatures..



Research Contributors and Funding

- Mr Cameron Heaton
- Miss Katie Kennedy
- Dr. L. Ferguson (Anglia Ruskin University)
- Dr R. Bradshaw (Sheffield Hallam University)
- Dr E. Patel (Shimadzu)
- Dr Gino Groeneveld (University of Amsterdam)
- Dr. R. Wolstenholme (Sheffield Hallam University)
- Prof. M.R. Clench (Sheffield Hallam University)
- Dr S. Bley (now LSBU), Mr Vaughn Sears (ex CAST, Home Office)**
- Dr H. Bandey, Mr Mark Seale, Dr R. McColm (Dstl)**
- Mr Neil Denison (Yorkshire and Humber RSSS)**
- Prof. M. de Puit (NFI, The Netherlands)**
- Dr G. Langenburg (Elite Forensic Services, US)**
- Michal Levin and Dr Ravell Bengiat (Israel Police)**
- Prof. Marcel dePuit (NFI)**
- Dr Charlie Bury, Catapult Medicine Discovery**

